



# Material Data Sheet

## General Information

### Headquarters

**TELIX Pharmaceuticals Ltd**  
55 Flemington Road, Suite 401  
North Melbourne, VIC 3051  
Australia

### Manufacturing Facility

**TELIX Innovations SA**  
Rue de Hermée, 255  
4040 Herstal  
Belgium

## Material Ordering Description and Overview

### For EZAG Generator

Optimized for 5 mL eluate of 0.1 M HCl



### Box of 5 Sterile Acetate Buffer Glass Vials, Config. A

ITEM CODE	PRODUCT REFERENCE	SHELF LIFE
900018	BOX-GSV-BF-AC-A	36 months

### For IRE Generator

Optimized for 1.1 mL eluate of 0.1 M HCl



### Box of 5 Sterile Acetate Buffer Glass Vials, Config. B

ITEM CODE	PRODUCT REFERENCE	SHELF LIFE
900019	BOX-GSV-BF-AC-B	36 months

### Sterile Vacuumed Glass Vial



### Box of 5 Sterile Vacuumed Glass Vials Compatible with both generators

ITEM CODE	PRODUCT REFERENCE	SHELF LIFE
900017	BOX-SVV-10G	24 months

### Pack of Accessories



### For EZAG Generator

- 5 Syringes of 10 mL
- 5 Swabable Extension Lines
- 15 Vented Swabable Vial Adaptors

ITEM CODE	SHELF LIFE
900020	24 months



### For IRE Generator

- 5 Syringes of 10 mL
- 5 Swabable Extension Lines
- 15 Non-Vented Swabable Vial Adaptors

ITEM CODE	SHELF LIFE
900021	24 months

Storage conditions for all materials above: *Do not freeze.*